

# Lymph-Pro

## Digital

by



*Innovating for Life*


# The Power of Compression Physiotherapy

Compression therapy is a cornerstone of modern physical rehabilitation, designed to stimulate circulation, reduce swelling, and promote the body's natural healing processes. By applying controlled, rhythmic pressure to the limbs, it effectively manages lymphedema, reduces venous insufficiency, and accelerates recovery post-surgery or injury. This non-invasive technique is essential for enhancing patient comfort, improving mobility, and restoring quality of life. Until now, accessing this level of therapeutic care required complex equipment and specialized knowledge. Health Technology is changing that narrative.

## Introducing **LYMPH-PRO** Digital

### Where Innovation Meets Simplicity

Welcome to the future of compression therapy. The Lymph Pro Digital, manufactured with precision in Egypt by Health Technology, is engineered to set a new standard for patient care and clinical efficiency. This device stands alone in the market with its integrated digital display and intelligently pre-programmed treatment protocols. We've transformed advanced medical technology into an exceptionally user-friendly experience, ensuring that both practitioners and patients can achieve optimal results with confidence and ease. It's not just a device; it's your partner in proactive health.



### Principle:

#### Mode A:

Typically common relaxing massage, effective hand massage type.

#### Mode B:

Lymphatic Edema or DVT prevention and treatment.

#### Mode C:

Mullti-relax deep massage and recovery for muscles.

#### Mode D:

Leg/ Arm/ Waist slimming

#### Mode E:

B Blood circulation and lymphatic drainage



## Key Features & Benefits

### 1. Intelligent Digital Display

**Benefit:** Complete control and transparency. The clear, backlit LCD screen provides real-time feedback on pressure settings, treatment mode, and session time, eliminating guesswork and ensuring the therapy is delivered precisely as prescribed.

### 2. Pre-Programmed Treatment Modes

**Benefit:** Simplicity and consistency. With a selection of expert-developed programs for different conditions (e.g., lymphedema, post-operative edema, sports recovery), users can initiate a professional-grade treatment at the touch of a button. This reduces setup time and ensures therapeutic integrity for every session.

### 3. Unmatched Reliability & Support

**Benefit:** Peace of mind. As with all Health Technology products, the Lymph Pro Digital is backed by our unwavering commitment to support. This includes:

**Guaranteed Maintenance:** Comprehensive service support to keep your device operational.

**Professional Training:** We ensure you or your staff know how to use the device effectively and safely.

**High-Quality Components:** Built to last in a clinical or home-care environment.

**Constant Spare Parts Availability:** Minimizing downtime and protecting your investment for the long term.

### 4. Superior Engineering & Comfort

**Benefit:** Effective and comfortable treatments. The device features precisely calibrated pressure control and comes with a range of high-quality, durable sleeves and attachments designed for maximum patient comfort and optimal therapeutic outcomes.

### 5. Exceptional Value

**Benefit:** We deliver this entire package of innovation, quality, and support at a very competitive price, making advanced compression therapy more accessible than ever before without compromising on performance.

## Applications

### 1. Prevention

- A. Lymphatic Edema
- B. Limb Paralysis
- C. Limb Convulsion

### 1. Improvement

- A. Fat dissolution
- B. Rheumatoid Arthritis
- C. DVT
- D. Varicocele Therapy

### 1. management

- A. Pregnant Woman Limb Management



+2 010 23448873



Health-Technology-Egypt



/HealthTechnologyCo



/co.healthtechnology



/HealthTechnology-medsup



## Technical Specifications

<b>Manufacturer:</b>	Health Technology, Egypt
<b>Model:</b>	Lymph Pro Digital
<b>Display:</b>	Digital LCD Screen
<b>Operation Modes:</b>	Multiple Pre-Programmed & Customizable Settings
<b>Pressure Range:</b>	30 - 220 mmHg
<b>Power Supply:</b>	220 V - 50 Hz
<b>Accessories:</b>	Left Leg Sleeve + Hose + Extension Right Leg Sleeve + Hose + Extension Arm Sleeve + Hose Abdomin Sleeve + Hose Dual Channel Air Flow Splitter

## The Health Technology Promise

When you choose the Lymph Pro Digital, you are investing in more than a device. You are partnering with a leader in medical equipment manufacturing in the Middle East. Our dedication to innovation, quality, and comprehensive customer care is built into every product we create. We stand behind our technology, so you can stand behind your treatment.

Ready to enhance your practice with  
Contact Us Now

**LYMPH-PRO**  
*Digital*

[www.health-technology.net](http://www.health-technology.net)

Main Office: Address: 18 El Shekh Ali Mahmoud – Mahkama Sqr. – Heliopolis, Cairo Egypt - Tel.: +202 264 341 80  
Factory: No. 317 Green Land Complex – Obour City, Egypt - Tel.: +202 448 127 08  
Email: [health@health-technology.net](mailto:health@health-technology.net)

